


Endodontic / Micro Suction Aspirators



Flexible Tip
(except 01BA)




overall length 6.0 inches

overall length 6.5 inches



- 0.5 mm
- 1.0 mm
- 1.7 mm
- 2.0 mm
- 2.3 mm



overall length 6.25 inches


overall length 6.25 inches




E04-90

E04

overall length 3.0 inches



04AT - 45° bend (shown) 04-90AT - 90° bend
overall length 8.5 inches



04T - 45° bend 04-90T - 90° bend (shown)
overall length 7.5 inches

overall length 5.5 inches



Disposable rubber tips
for the EY ENDO A
ERT (36 tips) **\$24.00**
RT (144 tips) **\$56.00**

Micro suction for root canal and root tip suctioning. Button provides finger rest and precise control of aspirator. The 02 and 03 feature bendable tubing for more accurate placement. Includes one cleaning stylet. Recommended cleaning brush: Extra Small.

Opening	Fits Into	Product #	Price
0.5 mm	HVE Cut-Off Valve	01BA	\$72.00
1.0 mm	HVE Cut-Off Valve	02BA	\$72.00
2.0 mm	HVE Cut-Off Valve	03BA	\$72.00

Opening	Fits Into	Product #	Price
1.0 mm	Surgical Tubing	02BAT	\$72.00
2.0 mm	Surgical Tubing	03BAT	\$72.00

Replacement cleaning stylet for endodontic aspirators.

For use with	Product #	Price	For use with	Product #	Price
01BA	ST-01BA	\$5.00	02BAT	ST-02BAT	\$5.00
02BA	ST-02BA	\$5.00	03BAT	ST-03BAT	\$5.00
03BA	ST-03BA	\$5.00			

Endodontic aspirator made of 1/8" diameter stainless steel. Recommended cleaning brush size (Pg. 25): Small.

Opening	Fits Into	Product #	Price
1.7 mm	HVE Cut-Off Valve	04A	\$30.00
2.3 mm	HVE Cut-Off Valve	05A	\$30.00

Opening	Fits Into	Product #	Price
1.7 mm	HVE Cut-Off Valve	04-90A	\$30.00

Endodontic aspirator made of 1/8" diameter stainless steel. Recommended cleaning brush size (Pg. 25): Small.

Opening	Fits Into	Product #	Price
1.7 mm	HVE Cut-Off Valve	E04	\$30.00
1.7 mm	HVE Cut-Off Valve	E04-90	\$30.00

Endodontic aspirator made of 1/8" diameter stainless steel. 8-1/2" long. Recommended cleaning brush size (Pg. 25): Small.

Opening	Fits Into	Product #	Price
1.7 mm	Surgical Tubing	04AT	\$32.00
1.7 mm	Surgical Tubing	04-90AT	\$32.00

Same as above but 7-1/2" long. Recommended cleaning brush size (Pg. 25): Small.

Opening	Fits Into	Product #	Price
1.7 mm	Surgical Tubing	04T	\$32.00
1.7 mm	Surgical Tubing	04-90T	\$32.00

The EY ENDO A aspirator uses disposable rubber tips and features a unique screw-off head which allows for easy cleaning. May be used in a saliva ejector or high volume valve.

Recommended cleaning brush size (Pg. 25): Medium.

Fits Into	Product #	Price
Saliva Ejector or High Volume Valve	EY ENDO A	\$69.00