

# Suction Valves & Accessories



Plug for valve quick disconnect. Can be used to plug the FLIP-A, FLIP-D, or FLIP-SE disconnect while the valve is being cleaned.

**HVE-PLUG \$9.95**



Replacement Screw On Head for FLIP-SE.

**SSE-HEAD \$18.00**



Autoclavable steel filter for FLIP-SE.

**SEFB \$6.00**

Screw On Head

Disposable plastic filters for FLIP-SE. Bag of 100

**P-SEFB \$9.00**



Replacement O'Ring set for the Flip-A, Flip-D valves.

**ROS-F \$4.50**



Replacement O'Ring set for the Flip-SE valve.

**ROS-FSE \$4.50**



The autoclavable Flip-A cut-off valve features easy on/off flip action, quick disconnect from the suction tubing, and an easy to remove and clean flip valve. Fits standard 1/2" vacuum tubing. Includes cleaning brush.

Fits Into	Product #	Price
Standard 1/2" Vacuum Tubing	<b>FLIP-A</b>	<b>\$54.00</b>

The autoclavable Flip-D cut-off offers the same features as the Flip-A, but with a shorter overall length. Fits standard 1/2" vacuum tubing. Includes cleaning brush.

Fits Into	Product #	Price
Standard 1/2" Vacuum Tubing	<b>FLIP-D</b>	<b>\$54.00</b>

The flip style autoclavable saliva ejector valve features a quick disconnect and easy to remove rotary for easy cleaning. The autoclavable silicone head screws on and off for a more secure fit.

Fits Into	Product #	Price
Saliva Ejector Tubing	<b>FLIP-SE</b>	<b>\$45.00</b>

Surgical tubing quick disconnect available for the FLIP-A, FLIP-D, and FLIP-SE valves. Allows you to attach the valve to 1/4" surgical tubing.

Connects to Valves	Product #	Price
FLIP-A or FLIP-D	<b>FA-ST-DISC</b>	<b>\$10.00</b>
FLIP-SE	<b>SE-ST-DISC</b>	<b>\$8.00</b>

Valve Holders. The small valve holder is recommended for A-Type valves and saliva ejector valves. The neck portion of the valve hangs in the bracket as shown. The large valve holder is recommended for the FLIP-D valve. The valve is held by the quick disconnect as shown.

Holds Valves	Product #	Price
FLIP-A or FLIP-SE	<b>SVH</b>	<b>\$16.00</b>
FLIP-D	<b>L VH</b>	<b>\$16.00</b>

# Saliva Ejectors



overall length  
6.25 inches



High volume saliva ejector fits directly into a high volume cut-off valve. Four autoclavable silicone sleeves are provided with each saliva ejector. Additional silicone sleeves are available in packages of 12. Recommended cleaning brush size (Pg. 25): Medium.

**SEA \$30.00**

Autoclavable replacement sleeves. Package of 12 sleeves.

**RSS-SEA \$22.00**

This all stainless steel saliva ejector has a removable head for easy cleaning. Ample holes clear the mouth of saliva and helps prevent grabbing soft tissue. Recommended brush size (Pg. 25): Small.

**NSE \$39.00**



Replacement head for NSE.  
**RHNSE \$26.00**

overall length  
4.75 inches



# Big Graft Material Syringe

Autoclavable stainless steel Big Graft Material Syringe with 3/16" diameter tubing for sinus augmentation and tunneling.

**BGMS-316 \$89.00**



Autoclavable stainless steel Big Graft Material Syringe for sinus augmentation and tunneling.

**BGMS-STR \$89.00**



Autoclavable stainless steel Big Graft Material Syringe for sinus augmentation and tunneling.

**BGMS \$89.00**