

# Cleaning Brushes & Stylets

Specially designed for aspirator tips and saliva ejectors, these brushes come in various sizes to accommodate a wide range of tubular diameters. The white nylon bristles clean blood and debris from inside walls. The stainless steel wire prevents rusting. Made 12" long, these brushes will extend the full length of aspirator tips and are flexible enough to round 45 and 90 degree curvatures.



- XSCB** (1 dozen Extra Small Cleaning Brushes) **\$28.00 dz.**
- SCB** (1 dozen Small Cleaning Brushes) **\$28.00 dz.**
- MCB** (1 dozen Medium Cleaning Brushes) **\$28.00 dz.**
- LCB** (1 dozen Large Cleaning Brushes) **\$28.00 dz.**
- XLCB** (1 dozen Extra Large Cleaning Brushes) **\$28.00 dz.**
- SECB** (1 dozen SE Cleaning Brushes) **\$28.00 dz.**
- CB** (1 dozen Assorted Cleaning Brushes) **\$28.00 dz.**

Autoclavable Bur / Instrument brush with stainless steel and nylon bristles and an easy to grip plastic handle. Ideal for cleaning burs, aspirators and dental instruments.

**DH-BB \$14.95**

Autoclavable Bur / Instrument brush with stainless steel bristles and an easy to grip plastic handle. Ideal for cleaning burs, aspirators and dental instruments.

**BB \$9.95**

Cleaning stylets for surgical aspirators. Standard length is 7 1/2". Silver soldered tips on the ends of the stylets correspond to the orifice size of the aspirators. They are best used for quickly dislodging clogs during surgery.

- ST7-0** (1.0 mm Cleaning Stylet) **\$4.00**
- ST7-1** (1.5 mm Cleaning Stylet) **\$9.00**
- ST7-2** (2.5 mm Cleaning Stylet) **\$9.00**
- ST7-3** (3.0 mm Cleaning Stylet) **\$9.00**